

TECHNICAL DATA SHEET

PRODUCT FAMILY: Finger-lift

<u>COMPOSITION</u>	<u>Weight</u>	<u>Thickness</u>
Adhesive protected side	35 gsm	Total 80 micron
Tissue carrier	13 gsm	
Adhesive open side	35 gsm	
Liner	62gsm	55 micron

PRODUCT DESCRIPTION

A permanent tape with a tissue carrier. The adhesive is an aggressive synthetic rubber based system and the liner a two side siliconised release paper which overlaps at each edge to give a 'finger-lift' for easy liner removal.

APPLICATIONS

The original and industry standard finger-lift tissue tape has been used for continuous stationery, mailing and padded envelopes, studio/commercial artwork, photographic, Point of Sale and splicing. Especially useful for many products produced by print finishers.

ADHESIVE PROPERTIES

			<u>Test method</u>
Peel Adhesion to steel (24 hour dwell)	protected	1300 N/M	PSTC 101
	open	1300 N/M	PSTC 101
Shear Adhesion to steel (1 kg; 25mm x 25mm)	protected	100 hours	PSTC 107
	open	100 hours	PSTC 107

THERMAL PROPERTY GUIDELINES

Service temperature range (short term)	0° C to 70° C
Application temperature range	10° C to 50° C

FOOD CONTACT INFORMATION

Based on supplier information the components of the adhesive conform to the compositional requirements of FDA 175.105 for indirect food contact. Additional to this compliance there should also be a functional barrier between adhesive and the food. Such functional barrier and determination of suitability are the responsibility of the customer.

STORAGE

For guidance purposes we recommend that the product be used within 12 months of supply. The product should be stored in original closed cartons at non extreme ambient temperatures away from direct sources of heat and sunlight.

NOTE

Values given are average and should not be taken as a maximum or minimum for specification purposes. Specific technical details, application information & suggested uses available.

WARRANTY

All data relating to suitable uses and description information concerning our products are compiled from research and are believed to be reliable but are provided for guidance only. The Company is not bound contractually unless and until an order acknowledgement in writing is made and must be on the Company's standard Condition of Sales. It is recommended that the Purchasers independently determine the suitability of the products for their purposes. No salesman, representative or agent is authorised to give any guarantee, warranty or make any representation contrary to the foregoing